



XCOLLAR, XCOLLAR PLUS, NEXSPLINT AND NEXSPLINT PLUS SPECIFICATIONS

| | |
|---|---|
|  | |
| PRODUCT: | XCOLLAR |
| Presentation: | Flat configuration, individually packed in a printed thin polyethylene bag. 30 units per case |
| Unit Size: | L= 57.2cm W= 21.6cm H= 1.27cm L= 22.5 in W= 8.5in H= 0.5in |
| Unit Weight: | 227g (8 oz) |
| Case Dimensions: | L= 61.7cm W= 46.7cm H= 22.2cm L= 24.3in W= 18.4in H= 8.7in |
| Gross Weight: | 7.5 Kg. (16.5 lbs.) |

XCOLLAR PRODUCT ID: **A-0100**
COLOR: **CLEAR (TRANSPARENT)**

A-0120
BRIGHT ORANGE

| | |
|---|---|
|  | |
| <p>Each XCollar Plus is equipped with an Integrated Adjustable Head Restraint System, consisting of a set of Velcro straps to secure the patient to the backboard from lateral and vertical movement, and an Adjustable Occipital Support Pad that ensures proper patient alignment.</p> | |
| PRODUCT: | XCOLLAR PLUS |
| Presentation: | Flat configuration, individually packed in a printed thin polyethylene bag. 28 units per case |
| Unit Size: | L= 57.2cm W= 21.6cm H= 2.47cm L= 22.5 in W= 8.5in H=1.0in |
| Unit Weight: | 272g (9.6 oz) |
| Case Dimensions: | L= 74.5cm W= 55cm H= 21.5cm L= 29.5in W= 21.7in H= 8.5in |
| Gross Weight: | 8.5Kg. (18 lbs.) |

XCOLLAR PLUS PRODUCT ID: **A-H0100**
COLOR: **CLEAR (TRANSPARENT)**

A-H0120
BRIGHT ORANGE



NeXsplint is the compact, foldable version of the XCollar, packed in a heavy duty, vacuum sealed, polyethylene bag to protect the device in extreme conditions. It has been designed for vehicle extrication, tactical medicine, over the side, confined space, mountain and urban rescue.

NeXsplint can be stored in a space 40% smaller than the XCollar and conventional c-collars (length = 15.0 in).

| | | | |
|-------------------------|--|-------------------------------|-------------------|
| PRODUCT: | NEXSPLINT | | |
| Presentation: | Folded configuration, packed in thick polyethylene vacuum sealed package 32 units per case | | |
| Unit Size: | L= 38.1cm W= 22.9cm H= 3.8cm | L= 15.0 in | W= 9.0in H= 1.5in |
| Unit Weight: | 232g (8.2 oz) | | |
| Case Dimensions: | L= 46.5cm W= 42.0cm H= 39.5cm | L= 18.4in W= 16.6in H= 15.6in | |
| Gross Weight: | 8.10 Kg. (17.85 lbs.) | | |

NEXSPLINT PRODUCT ID: B-0100

B-0120

COLOR: CLEAR (TRANSPARENT)

BRIGHT ORANGE



NeXsplint Plus can be stored in a space 40% smaller than the XCollar and conventional c-collars (length = 15.0 in).

Each NeXsplint Plus equipped with an Integrated Adjustable Head Restraint System, consisting of a set of Velcro straps to secure the patient to the backboard from lateral and vertical movement, and an Adjustable Occipital Support Pad that ensures proper patient alignment.

| | | | |
|-------------------------|--|-------------------------------|-------------------|
| PRODUCT: | NEXSPLINT | | |
| Presentation: | Folded configuration, packed in thick polyethylene vacuum sealed package 22 units per case | | |
| Unit Size: | L= 38.1cm W= 22.9cm H= 5cm | L= 15.0 in | W= 9.0in H= 2.0in |
| Unit Weight: | 277g (9.8 oz) | | |
| Case Dimensions: | L= 46.5cm W= 42.0cm H= 39.5cm | L= 18.4in W= 16.6in H= 15.6in | |
| Gross Weight: | 7.05 Kg. (15.5 lbs.) | | |

NEXSPLINT PLUS PRODUCT ID: B-H0100

B-H0110

B-H0120

COLOR: CLEAR (TRANSPARENT)

MILITARY GREEN

BRIGHT ORANGE